

# SCANLAN® *LEGACY* Titanium Forceps



- *Tapered DeBakey Tips*

- *Platform Tips*

- *Platform Ring Tips*

- *Platform Tips with  
1 x 2 Micro-Teeth*

SCANLAN ★

INTERNATIONAL

*Surgical products designed and manufactured by the Scanlan family since 1921*

# SCANLAN® LEGACY Titanium Forceps

These new titanium forceps offer you a choice of 6 different tip styles to meet your surgical needs:

- Tapered DeBakey Tips
- Platform Ring Tips
- 0.5 mm Platform Tips
- 0.5 mm Platform Tips with 1 x 2 Micro-Teeth
- 1.0 mm Platform Tips
- 1.0 mm Platform Tips with 1 x 2 Micro-Teeth

The unique *Never Shear™* DUAL GUIDE™ design of the SCANLAN® LEGACY forceps, prevents the forceps tips from scissoring and shearing delicate tissue.

The elegant blue and silver colors readily identify these as SCANLAN® LEGACY titanium forceps.

The well-defined knurling and slightly larger diameter provide a very secure, comfortable, and tactile grasp.



## SCANLAN® LEGACY Titanium forceps

### 1.0 mm Tapered DeBakey Tips

3003-932	Straight tips, 7-1/8" / 18.1 cm
3003-934	Straight tips, 8-1/4" / 21 cm
* 3003-936	Straight tips, 9" / 23 cm
* 3003-938	Straight tips, 10" / 25 cm
* 3003-806	Bayonet, straight tips, 8-3/4" / 22 cm

### DIAMOND DUST™ Platform Tips

0.5 mm Tips	1.0 mm Tips	Description
* 3003-902	3003-912	Straight tips, 7-1/8" / 18.1 cm
3003-904	3003-914	Straight tips, 8-1/4" / 21 cm
3003-906	3003-916	Straight tips, 9" / 23 cm
* 3003-908	* 3003-918	Straight tips, 10" / 25 cm
* 3003-816	* 3003-818	Bayonet, straight tips, 8-3/4" / 22 cm

### 1.0 mm DIAMOND DUST™ Platform Ring Tips

3003-922	Straight tips, 7-1/8" / 18.1 cm
3003-924	Straight tips, 8-1/4" / 21 cm
* 3003-926	Straight tips, 9" / 23 cm
* 3003-928	Straight tips, 10" / 25 cm
* 3003-804	Bayonet, straight tips, 8-3/4" / 22 cm

### DIAMOND DUST™ Platform Tips with 1 x 2 Micro Teeth

0.5 mm Tips	1.0 mm Tips	Description
3003-770		Puskas™, straight tips, 7 1/8"/18.1cm
3003-774	3003-778	Puskas™, straight tips, 8-1/4" / 21 cm
3003-775		Puskas™, curved tips, 8-1/4" / 21 cm
3003-776	3003-779	Puskas™, angled tips, 8-1/4" / 21 cm
* 3003-780	* 3003-784	Straight tips, 9" / 23 cm
* 3003-782	* 3003-786	Straight tips, 10" / 25 cm
3003-810		Bayonet, straight tips, 8-3/4" / 22 cm

\* Available Later 2008

The *Never Shear™* DUAL GUIDE™ design was developed in cooperation with John D. Puskas, M.D.

Patent Pending on the *Never Shear™* DUAL GUIDE™

For more information or to place an order, please contact your local Scanlan representative, or call...

U.S. & Canada: 800-328-9458

International: 651-298-0997

FAX: 651-298-0018 • [www.scanlaninternational.com](http://www.scanlaninternational.com)



INTERNATIONAL

*Surgical products designed and manufactured by the Scanlan family since 1921*

ONE SCANLAN PLAZA • SAINT PAUL, MINNESOTA 55107 • U.S.A.

800-328-9458 • 651-298-0997 • FAX 651-298-0018 • [www.scanlaninternational.com](http://www.scanlaninternational.com)

MK-B-NSDG Rev. B

© 2008 Scanlan International, Inc.